

## True3 Digital Benefits

### True Offloading through Sequential Whole-Body Pressure Relief

“Active” vs “Reactive” Technology - Alternating pressure relief regardless of patient’s movements or position.  
Whole Body Sequential Pressure Reduction to ZERO every 7.5 minutes

### TRI-PHASIC Alternation to promote Reactive Hyperemia

Sequential, “rolling” of 3 different pressures promotes capillary micro-bursts to help maximize blood flow.

### Self-Contained Design – No Powerful Fan Blowing Air Constantly

Extremely Quiet for patient comfort and to support sleep - No added Heat to the room for patient comfort

### Infection Control

No Bacteria/Virus spread from air sucked into pump or airflow out of cells, across patient, and into room

### Fail Safe “3rd Row” of air cells Supports Patient even through Power Failures

Bottom rows of cells dedicated to patient protection via a dedicated fill tube with a one-way air flow valve.  
Air cannot escape even if power is lost or pump is disconnected or malfunctioning.

### 3-Cell Support for More Comfort, Less Pain

4 inflated cells surround each deflated cell so that the patient is constantly supported.  
Patient lies “On” the bed instead of being immersed “In” the bed for a more natural body alignment.

### Staff Convenience

No Special Linens or Expensive Breathable Incontinence Pads Needed.  
Simple Design – Strap Down Mattress, Attach Pump, Plug In & Turn On.  
Prominent CPR Release at handle attachment to Pump

### Reliability

FULL 1 YEAR WARRANTY against manufacturer defects

[www.truewoundcare.com](http://www.truewoundcare.com)

P.O. Box 780241 San Antonio, TX 78278 | Phone: 888.880.TRUE | Email: [info@truewoundcare.com](mailto:info@truewoundcare.com)